



PF -Xpat

Malaria Pf (HRP2/pLDH) Test



The Malaria Pf-Xpat kit (HRP2/pLDH) is a lateral flow immunochromatographic test, intended for the qualitative detection of *P. falciparum* specific histidine rich protein-2 (P.f HRP2) and *P. falciparum* specific pLDH (P.f pLDH) in human whole blood sample.

- Less Whole Blood. Only 5µl of whole blood from the fingertip
- Easy-to-use
- Test in 3 steps, receive results in 15 mins.

Technical Specifications

- Test Type Qualitative Rapid Diagnostic Test
- Package 25 or 40 tests/box
- Specimen Type Whole Blood
- Specimen Volume 5 µl
- Result Time 15 - 30 mins
- Shelf Life 24 months

Cat. No.	Product	Storage temperature	Test/Kit
M002-C0001 (25 Tests) M002-C0002(40Tests)	Pf Xpat Malaria P.f (HRP2/pLDH) Rapid Test	1 - 40° C	25T/kit 40T/kit

Microhaem Scientifics

Email: info@microhaem.co.ug
sales@microhaem.co.ug
Tel: +256-393-246401 / +256-313-306900
Plot 16 A-C, Martyrs Way Ntinda, Uganda

www.microhaem.co.ug



@Microhaem Scientifics

Doc. No: 04-2024FL01-V01



Issue Date: 22/04/2024



**“Innovation,
Localised”**

Pioneering the future of IVD manufacturing in Africa.
Towards a **HIV-free future**



Kwiq-Test

HIV-1/2 Whole Blood/Serum/Plasma Test

The Kwiq-Test HIV 1/2 is a lateral flow immunochromatographic test, intended for the qualitative detection of HIV 1/2 antibodies in human whole blood, serum or plasma sample.

- Less sample. Only 10 µl of whole blood, serum or plasma
- Easy-to-use
- Test in 3 steps, receive results in 15 mins.

Technical Specifications

- Test Type: Qualitative Rapid Diagnostic Test
- Package: 25 or 40 tests/box
- Specimen Type: Whole Blood
- Specimen Volume: 10 µl
- Result Time: 15 - 30 mins
- Shelf Life: 24 months

Cat. No.
M001-C0001 (25 Tests)
M001-C0002(40Tests)

Product
HIV-1/2 Whole Blood/Serum/Plasma Test

Storage temperature
2 - 30 °C

Test/Kit
25T/kit 40T/kit

Microhaem Scientifics

Email: info@microhaem.co.ug
sales@microhaem.co.ug
 Tel: +256-393-246401 / +256-313-306900
 Plot 16 A-C, Martyrs Way Ntinda, Uganda

www.microhaem.co.ug



@Microhaem Scientifics

Doc. No: 04-2024FL01-V01



Issue Date: 22/04/2024